

RGF1A Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP13696
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human protein . at AA range: 130-210
Mol wt	N/A
Species reactivity	Human, Rat, Mouse
Clonality	Polyclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	Ras-GEF domain-containing family member 1A
Synonyms	Ras-GEF domain-containing family member 1A

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

Guanine nucleotide exchange factor (GEF) for KRAS, HRAS, and NRAS (in vitro). Plays a role in cell migration.,Contains 1 N-terminal Ras-GEF domain.,Contains 1 Ras-GEF domain.,tissue specificity:Detected in brain and spinal cord. Highly expressed in a number of intrahepatic cholangiocarcinoma tissue biopsies.,

Recommended Dilution

WB: 1: 500 - 1: 2000

ELISA: 1: 5000 - 1: 20000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com