

RGL2 Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP13700
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human protein . at AA range: 130-210
Mol wt	N/A
Species reactivity	Human, Mouse
Clonality	Polyclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	Ral guanine nucleotide dissociation stimulator-like 2
Synonyms	Ral guanine nucleotide dissociation stimulator-like 2 ;RalGDS-like 2;RalGDS-like factor;Ras-associated protein RAB2L

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

Probable guanine nucleotide exchange factor. Putative effector of Ras and/or Rap. Associates with the GTP-bound form of Rap 1A and H-Ras in vitro.,Contains 1 N-terminal Ras-GEF domain.,Contains 1 Ras-associating domain.,Contains 1 Ras-GEF domain.,

Recommended Dilution

WB: 1: 500 - 1: 2000

ELISA: 1: 5000 - 1: 20000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com