

S27A6 Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP14008
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human protein . at AA range: 270-350
Mol wt	N/A
Species reactivity	Human, Rat, Mouse
Clonality	Polyclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	Long-chain fatty acid transport protein 6
Synonyms	Long-chain fatty acid transport protein 6 ;FATP-6;Fatty acid transport protein 6;Fatty-acid-coenzyme A ligase, very long-chain 2;Solute carrier family 27 member 6;Very long-chain acyl-CoA synt

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

This gene encodes a member of the fatty acid transport protein family (FATP). FATPs are involved in the uptake of long-chain fatty acids and have unique expression patterns. Alternatively spliced transcript variants encoding the same protein have been found for this gene.

Recommended Dilution

WB: 1: 500 - 1: 2000

ELISA: 1: 5000 - 1: 20000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com