

S6A13 Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP14021
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from part region of human protein
Mol wt	N/A
Species reactivity	Human, Rat, Mouse
Clonality	Polyclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	Sodium- and chloride-dependent GABA transporter 2
Synonyms	Sodium- and chloride-dependent GABA transporter 2 ;GAT-2;Solute carrier family 6 member 13

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

Terminates the action of GABA by its high affinity sodium-dependent reuptake into presynaptic terminals.,Belongs to the sodium:neurotransmitter symporter (SNF) family.,tissue specificity:Expressed in brain, kidney, lung, liver and testis.,

Recommended Dilution

WB: 1: 500 - 1: 2000

ELISA: 1: 5000 - 1: 20000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com