

SCAFB Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP14059
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from part region of human protein
Mol wt	N/A
Species reactivity	Human, Rat, Mouse
Clonality	Polyclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	Protein SCAF11
Synonyms	Protein SCAF11 ;CTD-associated SR protein 11;Renal carcinoma antigen NY-REN-40;SC35-interacting protein 1;SR-related and CTD-associated factor 11;SRSF2-interacting protein;Serine/arginine-ri

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

Plays a role in pre-mRNA alternative splicing by regulating spliceosome assembly.,PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,sequence Contaminating sequence. Potential poly-A sequence.,subunit:Interacts with SFRS2, U2AF2 and SNRNP70.,tissue specificity:Widely expressed.,

Recommended Dilution

WB: 1: 500 - 1: 2000

ELISA: 1: 5000 - 1: 20000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com