

Secretogranin-1 Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP14103
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human Secretogranin-1
Mol wt	N/A
Species reactivity	Human, Mouse, Rat
Clonality	Polyclonal
Recommended application	IHC-p, IF, WB
Concentration	1 mg/ml
Full name	Secretogranin-1
Synonyms	Secretogranin-1; Secretogranin-1; Chromogranin-B; CgB; Secretogranin I; SgI; GAWK peptide; CCB peptide;

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

This gene encodes a tyrosine-sulfated secretory protein abundant in peptidergic endocrine cells and neurons. This protein may serve as a precursor for regulatory peptides.

Recommended Dilution

WB: 1: 500 - 1: 2000

IHC-p: 1: 50 - 1: 200

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com