

SEM4B Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP14107
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from part region of human protein
Mol wt	N/A
Species reactivity	Human, Rat, Mouse
Clonality	Polyclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	Semaphorin-4B
Synonyms	Semaphorin-4B

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

Inhibits axonal extension by providing local signals to specify territories inaccessible for growing axons.,Belongs to the semaphorin family.,Contains 1 Ig-like C2-type (immunoglobulin-like) domain.,Contains 1 PSI domain.,Contains 1 Sema domain.,

Recommended Dilution

WB: 1: 500 - 1: 2000

ELISA: 1: 5000 - 1: 20000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com