

ALPK3 Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP14195
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human protein . at AA range: 730-810
Mol wt	N/A
Species reactivity	Human, Mouse
Clonality	Polyclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	Alpha-protein kinase 3
Synonyms	Alpha-protein kinase 3 ;EC 2.7.11.-;Muscle alpha-protein kinase

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

Kinases that recognize phosphorylation sites in which the surrounding peptides have an alpha-helical conformation. Plays a role in cardiomyocyte differentiation.,Belongs to the protein kinase superfamily. Alpha-type protein kinase family. ALPK subfamily.,Contains 1 alpha-type protein kinase domain.,Contains 1 Ig-like (immunoglobulin-like) domain.,

Recommended Dilution

WB: 1: 500 - 1: 2000

ELISA: 1: 5000 - 1: 20000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com