

ATG13 Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP14386
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Recombinant Protein of ATG13
Mol wt	N/A
Species reactivity	Human, Mouse, Rat
Clonality	Polyclonal
Recommended application	IHC-p, IF
Concentration	1 mg/ml
Full name	Autophagy-related protein 13
Synonyms	Autophagy-related protein 13; Autophagy-related protein 13

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

The protein encoded by this gene is an autophagy factor and a target of the TOR kinase signaling pathway. The encoded protein is essential for autophagosome formation and mitophagy.

Recommended Dilution

IHC-p: 1: 50 - 1: 300

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com