

SPRL1 Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP14434
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from part region of human protein
Mol wt	N/A
Species reactivity	Human, Rat, Mouse
Clonality	Polyclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	SPARC-like protein 1
Synonyms	SPARC-like protein 1 ;High endothelial venule protein;Hevin;MAST 9

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

Belongs to the SPARC family.,Contains 1 EF-hand domain.,Contains 1 follistatin-like domain.,Contains 1 Kazal-like domain.,tissue specificity:Highly expressed in lymph node, brain, heart, lung, skeletal muscle, ovary, small intestine, and colon, with lower levels in placenta, pancreas, testis, spleen, and thymus, and no expression in kidney, liver, and peripheral blood leukocytes.,

Recommended Dilution

WB: 1: 500 - 1: 2000

ELISA: 1: 5000 - 1: 20000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com