

TBA8 Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP14783
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from part region of human protein
Mol wt	N/A
Species reactivity	Human, Mouse, Rat
Clonality	Polyclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	Tubulin alpha-8 chain
Synonyms	Tubulin alpha-8 chain ;Alpha-tubulin 8;Tubulin alpha chain-like 2

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

TUBA8 encodes a member of the alpha tubulin protein family. Alpha tubulins are one of two core protein families (alpha and beta tubulins) that heterodimerize and assemble to form microtubules. Mutations in this gene are associated with polymicrogyria and optic nerve hypoplasia. Alternate splicing results in multiple transcript variants. TUBA8 (Tubulin Alpha 8) is a Protein Coding gene. Diseases associated with TUBA8 include Cortical Dysplasia, Complex, With Other Brain Malformations 8 and Optic Nerve Hypoplasia. Among its related pathways are Development Slit-Robo signaling and Cytoskeleton remodeling Neurofilaments.

Recommended Dilution

WB: 1: 500 - 1: 2000

ELISA: 1: 5000 - 1: 20000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com