

TBD Polyclonal Antibody

Description

| | |
|--------------------------------|---|
| Product type | Primary Antibody |
| Code | BT-AP14789 |
| Host | Rabbit |
| Isotype | IgG |
| Size | 100ul, 50ul, 20ul |
| Immunogen | Synthesized peptide derived from part region of human protein |
| Mol wt | N/A |
| Species reactivity | Human, Rat, Mouse |
| Clonality | Polyclonal |
| Recommended application | WB, ELISA |
| Concentration | 1 mg/ml |
| Full name | Tubulin delta chain |
| Synonyms | Tubulin delta chain ;Delta-tubulin |

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

alternative products:Experimental confirmation may be lacking for some isoforms,In the elongating spermatid it is associated with the manchette, a specialized microtubule system present during reshaping of the sperm head.,Belongs to the tubulin family.,subcellular location:Associated with centrioles. Both cytoplasmic and nuclear.,

Recommended Dilution

WB: 1: 500 - 1: 2000

ELISA: 1: 5000 - 1: 20000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com