

TBD2A Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP14790
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from part region of human protein
Mol wt	N/A
Species reactivity	Human, Rat, Mouse
Clonality	Polyclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	TBC1 domain family member 2A
Synonyms	TBC1 domain family member 2A ;Armus;Prostate antigen recognized and identified by SEREX 1;PARIS-1

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

May act as a GTPase-activating protein for Rab family protein(s).,Contains 1 PH domain.,Contains 1 Rab-GAP TBC domain.,tissue specificity:Expressed in a broad range of tissues. Differentially expressed in prostate normal and cancer cells.,

Recommended Dilution

WB: 1: 500 - 1: 2000

ELISA: 1: 5000 - 1: 20000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com