

TM9S1 Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP15018
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from part region of human protein
Mol wt	N/A
Species reactivity	Human, Mouse, Rat
Clonality	Polyclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	Transmembrane 9 superfamily member 1
Synonyms	Transmembrane 9 superfamily member 1 ;MP70 protein family member;hMP70

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

Belongs to the nonspanin (TM9SF) family.,tissue specificity:Expressed in lung, pancreas, kidney, liver, placenta, skeletal muscle, heart and brain. The amount in skeletal muscle, heart and brain were considerably lower than in the other tissues.,

Recommended Dilution

WB: 1: 500 - 1: 2000

ELISA: 1: 5000 - 1: 20000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com