

TPTE2 Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP15102
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human protein . at AA range: 300-380
Mol wt	N/A
Species reactivity	Human, Rat, Mouse
Clonality	Polyclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	Phosphatidylinositol 3,4,5-trisphosphate 3-phosphatase TPTE2
Synonyms	Phosphatidylinositol 3,4,5-trisphosphate 3-phosphatase TPTE2 ;EC 3.1.3.67;Lipid phosphatase TPIP;TPTE and PTEN homologous inositol lipid phosphatase

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

TPIP is a member of a large class of membrane-associated phosphatases with substrate specificity for the 3-position phosphate of inositol phospholipids.

Recommended Dilution

WB: 1: 500 - 1: 2000

ELISA: 1: 5000 - 1: 20000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com