

CTSRB Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP15121
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from part region of human protein
Mol wt	N/A
Species reactivity	Human, Rat, Mouse
Clonality	Polyclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	Cation channel sperm-associated protein subunit beta
Synonyms	Cation channel sperm-associated protein subunit beta ;CatSper-beta

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

Probably involved in sperm cell hyperactivation via its association with CATSPER1. Sperm cell hyperactivation is needed for sperm motility which is essential late in the preparation of sperm for fertilization.,subunit:Component of a complex, at least composed of CATSPER1, HSPA1B and CATSPERB.,

Recommended Dilution

WB: 1: 500 - 1: 2000

ELISA: 1: 5000 - 1: 20000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com