

TRIM72 Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP15148
Host	Rabbit
Isotype	IgG
Size	20ul, 50ul, 100ul
Immunogen	Recombinant Protein of TRIM72
Mol wt	N/A
Species reactivity	Human, Mouse, Rat
Clonality	Polyclonal
Recommended application	WB, IHC-p, IF
Concentration	N/A
Full name	Tripartite motif-containing protein 72
Synonyms	Tripartite motif-containing protein 72; Tripartite motif-containing protein 72; Mitsugumin-53; Mg53

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

Belongs to the TRIM/RBCC family.,Contains 1 B box-type zinc finger.,Contains 1 B30.2/SPRY domain.,Contains 1 RING-type zinc finger.,


Recommended Dilution

WB: 1: 1000


IHC: 1: 200 - 1: 500

Not yet tested in other applications.

Images



Immunohistochemical analysis of paraffin-embedded Mouse Pancreas Tissue using TRIM72 Polyclonal Antibody.



Western blot analysis of Mouse skeletal muscle using TRIM72 Polyclonal Antibody. Secondary antibody was diluted at 1:20000

Storage

-20°C for 1 year