

TSH3 Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP15209
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from part region of human protein
Mol wt	N/A
Species reactivity	Human, Rat, Mouse
Clonality	Polyclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	Teashirt homolog 3
Synonyms	Teashirt homolog 3 ;Zinc finger protein 537

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

This gene encodes a zinc-finger transcription factor that regulates smooth muscle cell differentiation in the developing urinary tract. Consistent with this role, mice in which this gene has been inactivated exhibit abnormal gene expression in urinary tract smooth muscle cell precursors and kidney defects including hydronephrosis. The encoded transcription factor comprises a gene silencing complex that inhibits caspase expression. Reduced expression of this gene and consequent caspase upregulation may be correlated with progression of Alzheimer's disease in human patients.

Recommended Dilution

WB: 1: 500 - 1: 2000

ELISA: 1: 5000 - 1: 20000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year