

UBP37 Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP15325
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from part region of human protein
Mol wt	N/A
Species reactivity	Human, Rat, Mouse
Clonality	Polyclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	Ubiquitin carboxyl-terminal hydrolase 37
Synonyms	Ubiquitin carboxyl-terminal hydrolase 37 ;EC 3.4.19.12;Deubiquitinating enzyme 37;Ubiquitin thioesterase 37;Ubiquitin-specific-processing protease 37

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

Catalytic activity:Ubiquitin C-terminal thioester + H(2)O = ubiquitin + a thiol.,Belongs to the peptidase C19 family.,Contains 3 UIM (ubiquitin-interacting motif) repeats.,tissue specificity:Expressed in brain and prostate.,

Recommended Dilution

WB: 1: 500 - 1: 2000

ELISA: 1: 5000 - 1: 20000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com