

UBP54 Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP15329
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from part region of human protein
Mol wt	N/A
Species reactivity	Human, Rat, Mouse
Clonality	Polyclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	Inactive ubiquitin carboxyl-terminal hydrolase 54
Synonyms	Inactive ubiquitin carboxyl-terminal hydrolase 54 ;Inactive ubiquitin-specific peptidase 54

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

Although the active site residues are conserved, it lacks the conserved His residue which is normally found 9 residues before the catalytic His. Has no peptidase activity. Belongs to the peptidase C19 family. tissue specificity: Weakly expressed in a few tissues.

Recommended Dilution

WB: 1: 500 - 1: 2000

ELISA: 1: 5000 - 1: 20000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com