

ECFP-Tag Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP15360
Host	Mouse
Isotype	IgG
Size	20ul, 50uL, 100ul
Immunogen	Recombinant Protein of ECFP-Tag
Mol wt	N/A
Species reactivity	Species independent
Clonality	Polyclonal
Recommended application	WB

Concentration

Full name

Synonyms

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background


Cyan Fluorescent Protein (CFP) is a genetic mutant of green fluorescent protein (GFP) originally derived from the jellyfish *Aequorea victoria*.

Recommended Dilution

WB: 1: 2000 - 5000

Not yet tested in other applications.

Images



Western blot analysis of ECFP Recombinant protein, diluted at 1:5000; Secondary antibody was diluted at 1:20000

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com