

## VASN Polyclonal Antibody

### Description

|                                |   |
|--------------------------------|---|
| <b>Product type</b>            | Primary Antibody  |
| <b>Code</b>                    | BT-AP15414  |
| <b>Host</b>                    | Rabbit  |
| <b>Isotype</b>                 | IgG   |
| <b>Size</b>                    | 100ul, 50ul, 20ul   |
| <b>Immunogen</b>               | Synthesized peptide derived from part region of human protein |
| <b>Mol wt</b>                  | N/A   |
| <b>Species reactivity</b>      | Human, Mouse  |
| <b>Clonality</b>               | Polyclonal  |
| <b>Recommended application</b> | WB, ELISA   |
| <b>Concentration</b>           | 1 mg/ml   |
| <b>Full name</b>               | Vasorin   |
| <b>Synonyms</b>                | Vasorin ;Protein slit-like 2                                  |

**This product is for research use only, not for use in human, therapeutic or diagnostic procedure.**

### Background

May act as an inhibitor of TGF-beta signaling.,PTM:N-glycosylated.,Contains 1 EGF-like domain.,Contains 1 fibronectin type-III domain.,Contains 11 LRR (leucine-rich) repeats.,subunit:Interacts with TGFB1, TGFB2 and TGFB3.,tissue specificity:Expressed at highest levels in aorta, at intermediate levels in kidney and placenta and at lowest levels in brain, heart, liver, lung and skeletal muscle. Within the aorta, the strongest expression is found in the tunica media of the proximal ascending aorta, the descending thoracic aorta, the abdominal aorta and the coronary arteries. Within the kidney, expression is found in the interstitial cells.,

### Recommended Dilution

WB: 1: 500 - 1: 2000

ELISA: 1: 5000 - 1: 20000

Not yet tested in other applications.

### Images

No images.

### Storage

-20°C for 1 year