

F176A Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP15496
Host	Rabbit
Isotype	IgG
Size	20ul, 50ul, 100ul
Immunogen	Synthesized peptide derived from human protein . at AA range: 91-140
Mol wt	N/A
Species reactivity	Human, Mouse
Clonality	Polyclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	Protein FAM176A
Synonyms	Protein FAM176A ;Transmembrane protein 166

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

Acts as a regulator of programmed cell death, mediating both autophagy and apoptosis.,Belongs to the FAM176 family.,tissue specificity:Expressed in lung, kidney, liver, pancreas, placenta, but not in heart and skeletal muscle.,


Recommended Dilution

WB: 1: 500 - 1: 2000

ELISA: 1: 5000 - 1: 20000

Not yet tested in other applications.

Images



Western blot analysis of lysates from KB cells, primary antibody was diluted at 1:1000, 4°C overnight

Storage

-20°C for 1 year