

ZBT17 Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP15609
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human protein . at AA range: 160-240
Mol wt	N/A
Species reactivity	Human, Mouse
Clonality	Polyclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	Zinc finger and BTB domain-containing protein 17
Synonyms	Zinc finger and BTB domain-containing protein 17 ;Myc-interacting zinc finger protein 1;Miz-1;Zinc finger protein 151;Zinc finger protein 60

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

This gene encodes a zinc finger protein involved in the regulation of c-myc. The symbol MIZ1 has also been associated with PIAS2 which is a different gene located on chromosome 18.

Recommended Dilution

WB: 1: 500 - 1: 2000

ELISA: 1: 5000 - 1: 20000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com