

ZFH2 Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP15627
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from part region of human protein
Mol wt	N/A
Species reactivity	Human, Rat, Mouse
Clonality	Polyclonal
Recommended application	IHC-p, IF
Concentration	1 mg/ml
Full name	Zinc finger homeobox protein 2
Synonyms	Zinc finger homeobox protein 2 ;Zinc finger homeodomain protein 2;ZFH-2

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

Recommended Dilution

IHC-p: 1: 50 - 1: 300

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year