

ZN238 Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP15650
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from part region of human protein
Mol wt	N/A
Species reactivity	Human, Rat, Mouse
Clonality	Polyclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	Zinc finger protein 238
Synonyms	Zinc finger protein 238 ;58 kDa repressor protein;Transcriptional repressor RP58;Translin-associated zinc finger protein 1;TAZ-1;Zinc finger and BTB domain-containing protein 18;Zinc finger

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

This gene encodes a C2H2-type zinc finger protein which acts a transcriptional repressor of genes involved in neuronal development. The encoded protein recognizes a specific sequence motif and recruits components of chromatin to target genes. Alternative splicing results in multiple transcript variants.

Recommended Dilution

WB: 1: 500 - 1: 2000

ELISA: 1: 5000 - 1: 20000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year