

## GOT1B Polyclonal Antibody

### Description

Product type	Primary Antibody
Code	BT-AP15758
Host	Rabbit
Isotype	IgG
Size	20ul, 50ul, 100ul
Immunogen	Synthesized peptide derived from human protein . at AA range: 70-150
Mol wt	N/A
Species reactivity	Human, Mouse
Clonality	Polyclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	Vesicle transport protein GOT1B
Synonyms	Vesicle transport protein GOT1B ;Golgi transport 1 homolog B;Putative NF-kappa-B-activating protein 470;hGOT1a

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

### Background


#### Recommended Dilution

WB: 1: 500 - 1: 2000

ELISA: 1: 5000 - 1: 20000

Not yet tested in other applications.

### Images



### Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: [save@bt-laboratory.com](mailto:save@bt-laboratory.com) | [www.bt-laboratory.com](http://www.bt-laboratory.com)