

GST-Tag Polyclonal Antibody

Description

| | |
|--------------------------------|--------------------------------|
| Product type | Primary Antibody |
| Code | BT-AP15800 |
| Host | Rabbit |
| Isotype | IgG |
| Size | 100ul, 50uL, 20ul |
| Immunogen | Recombinant Protein of GST-Tag |
| Mol wt | N/A |
| Species reactivity | Species independent |
| Clonality | Polyclonal |
| Recommended application | WB |

Concentration

Full name

Synonyms

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

GST (Glutathione S-Transferase) is a widely used protein tag encoded by the *Schistosoma japonicum*. GST provides both an easily detectable tag and a simple purification process with little effect on the biological function of the protein of interest. Antibodies to GST are useful for checking protein expression both in plaques and on Western blots as well as for immunoaffinity purification of proteins expressed in cells.

Recommended Dilution

WB: 1: 5000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com