

IKKBeta Monoclonal Antibody

Description

Product type	Primary Antibody
Code	BT-MCA0108
Host	Mouse
Isotype	IgG
Size	50ul, 100ul
Immunogen	Purified recombinant fragment of IKKBeta expressed in E. Coli.
Mol wt	N/A
Species reactivity	Human
Clonality	Monoclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	Inhibitor of nuclear factor kappa-B kinase subunit beta
Synonyms	IKKB; IKK β ; Inhibitor of nuclear factor kappa-B kinase subunit beta; I-kappa-B-kinase beta; IKK-B; IKK-beta; Ikbkb; I-kappa-B kinase 2; IKK2; Nuclear factor NF-kappa-B inhibitor kinase beta; NFKB1KB

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

The protein encoded by this gene phosphorylates the inhibitor in the inhibitor/NF-kappa-B complex, causing dissociation of the inhibitor and activation of NF-kappa-B. The encoded protein itself is found in a complex of proteins. Several transcript variants, some protein-coding and some not, have been found for this gene.


Recommended Dilution

ELISA: 1:10000

WB: 1:500 - 1:2000

Not yet tested in other applications.

Images



Western Blot analysis using IKKBeta Monoclonal antibody against truncated IKKBeta recombinant protein (1) and K562 cell lysate (2).

Storage

-20°C for one year