

## Blk Monoclonal Antibody

### Description

Product type	Primary Antibody
Code	BT-MCA0216
Host	Mouse
Isotype	IgG
Size	50ul, 100ul
Immunogen	Purified recombinant fragment of Blk expressed in E. Coli.
Mol wt	N/A
Species reactivity	Human
Clonality	Monoclonal
Recommended application	WB, IHC-p, IF, ELISA
Concentration	1 mg/ml
Full name	Tyrosine-protein kinase Blk
Synonyms	BLK; Tyrosine-protein kinase Blk; B lymphocyte kinase; p55-Blk

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

### Background

N/A

### Recommended Dilution


ELISA: 1:10000

IHC: 1:200 - 1:1000


WB: 1:500 - 1:2000

Not yet tested in other applications.

### Images



Western Blot analysis using Blk Monoclonal antibody against truncated BLK recombinant protein Raji cell lysate.



Immunohistochemistry analysis of paraffin-embedded human breast tissue (A), lymph tissue (B) and skin carcinoma (C), showing membrane localization with DAB staining using Blk Monoclonal antibody.

### Storage

-20°C for one year