

CD3 Monoclonal Antibody

Description

Product type	Primary Antibody
Code	BT-MCA0294
Host	Mouse
Isotype	IgG
Size	50ul, 100ul
Immunogen	Purified recombinant fragment of human CD3 expressed in E. Coli.
Mol wt	N/A
Species reactivity	Human
Clonality	Monoclonal
Recommended application	FCM, ELISA
Concentration	1 mg/ml
Full name	T-cell surface glycoprotein CD3 zeta chain
Synonyms	CD247; CD3Z; T3Z; TCRZ; T-cell surface glycoprotein CD3 zeta chain; T-cell receptor T3 zeta chain; CD antigen CD247

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

The protein encoded by this gene is T-cell receptor zeta, which together with T-cell receptor alpha/beta and gamma/delta heterodimers, and with CD3-gamma, -delta and -epsilon, forms the T-cell receptor-CD3 complex. The zeta chain plays an important role in coupling antigen recognition to several intracellular signal-transduction pathways. Low expression of the antigen results in impaired immune response. Two alternatively spliced transcript variants encoding distinct isoforms have been found for this gene.


Recommended Dilution

ELISA: 1:10000

FC: 1:200 - 1:400

Not yet tested in other applications.

Images



Flow cytometric analysis CD3 expressed in human peripheral blood cells using CD3 Monoclonal antibody.

Storage

-20°C for one year