

DBC-1 Monoclonal Antibody

Description

Product type	Primary Antibody
Code	BT-MCA0437
Host	Mouse
Isotype	IgG
Size	50ul, 100ul
Immunogen	Purified recombinant human DBC-1 (N-terminus) protein fragments expressed in E.coli.
Mol wt	N/A
Species reactivity	Human,Mouse,Rat,Bovine,Dog,Pig
Clonality	Monoclonal
Recommended application	WB
Concentration	1 mg/ml
Full name	Protein KIAA1967
Synonyms	KIAA1967; DBC1; DBIRD complex subunit KIAA1967; Deleted in breast cancer gene 1 protein; DBC-1; DBC.1; p30 DBC

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background


N/A

Recommended Dilution

WB: 1:1000 - 1:2000

Not yet tested in other applications.

Images



Western Blot analysis using DBC-1 Monoclonal antibody against HeLa cell lysate.

Storage

-20°C for one year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com