

## FoxP2 Monoclonal Antibody

### Description

|                                |  |
|--------------------------------|--|
| <b>Product type</b>            | Primary Antibody   |
| <b>Code</b>                    | BT-MCA0568   |
| <b>Host</b>                    | Mouse  |
| <b>Isotype</b>                 | IgG  |
| <b>Size</b>                    | 50ul, 100ul  |
| <b>Immunogen</b>               | Purified recombinant fragment of human FoxP2 expressed in E. Coli.   |
| <b>Mol wt</b>                  | N/A  |
| <b>Species reactivity</b>      | Human  |
| <b>Clonality</b>               | Monoclonal   |
| <b>Recommended application</b> | WB, ELISA  |
| <b>Concentration</b>           | 1 mg/ml  |
| <b>Full name</b>               | Forkhead box protein P2  |
| <b>Synonyms</b>                | FOXP2; CAGH44; TNRC10; Forkhead box protein P2; CAG repeat protein 44; Trinucleotide repeat-containing gene 10 protein |

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

### Background

This gene encodes a member of the forkhead/winged-helix (FOX) family of transcription factors. It is expressed in fetal and adult brain as well as in several other organs such as the lung and gut. The protein product contains a FOX DNA-binding domain and a large polyglutamine tract and is an evolutionarily conserved transcription factor, which may bind directly to approximately 300 to 400 gene promoters in the human genome to regulate the expression of a variety of genes. This gene is required for proper development of speech and language regions of the brain during embryogenesis, and may be involved in a variety of biological pathways and cascades that may ultimately influence language development. Mutations in this gene cause speech-language disorder 1 (SPCH1), also known as autosomal dominant speech and language disorder with orofacial dyspraxia. Multiple alternative transcripts encoding different isoforms


### Recommended Dilution

ELISA: 1:10000

WB: 1:500 - 1:2000

Not yet tested in other applications.

### Images



Western Blot analysis using FoxP2 Monoclonal antibody against recombinant protein.

### Storage

-20°C for one year