

Fyb Monoclonal Antibody

Description

Product type	Primary Antibody
Code	BT-MCA0573
Host	Mouse
Isotype	IgG
Size	50ul, 100ul
Immunogen	Purified recombinant human Fyb protein fragments expressed in E.coli.
Mol wt	N/A
Species reactivity	Human,Rat
Clonality	Monoclonal
Recommended application	WB
Concentration	1 mg/ml
Full name	FYN-binding protein
Synonyms	FYB; SLAP130; FYN-binding protein; Adhesion and degranulation promoting adaptor protein; ADAP; FYB-120; 130; p120; p130; FYN-T-binding protein; SLAP-130; SLP-76-associated phosphoprotein

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background


The protein encoded by this gene is an adapter for the FYN protein and LCP2 signaling cascades in T-cells. The encoded protein is involved in platelet activation and controls the expression of interleukin-2. Three transcript variants encoding different isoforms have been found for this gene.

Recommended Dilution

WB: 1:1000 - 1:2000

Not yet tested in other applications.

Images



Western Blot analysis using Fyb Monoclonal antibody against Jurkat, K562. TF-1 cell lysate.

Storage

-20°C for one year