

GTBP Monoclonal Antibody

Description

Product type	Primary Antibody
Code	BT-MCA0647
Host	Mouse
Isotype	IgG
Size	50ul, 100ul
Immunogen	Purified recombinant fragment of GTBP expressed in E. Coli.
Mol wt	N/A
Species reactivity	Human
Clonality	Monoclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	DNA mismatch repair protein Msh6
Synonyms	MSH6; GTBP; DNA mismatch repair protein Msh6; hMSH6; G; T mismatch-binding protein; GTBP; GTMBP; MutS-alpha 160 kDa subunit; p160

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

This gene encodes a member of the DNA mismatch repair MutS family. In E. coli, the MutS protein helps in the recognition of mismatched nucleotides prior to their repair. A highly conserved region of approximately 150 aa, called the Walker-A adenine nucleotide binding motif, exists in MutS homologs. The encoded protein heterodimerizes with MSH2 to form a mismatch recognition complex that functions as a bidirectional molecular switch that exchanges ADP and ATP as DNA mismatches are bound and dissociated. Mutations in this gene may be associated with hereditary nonpolyposis colon cancer, colorectal cancer, and endometrial cancer. Transcripts variants encoding different isoforms have been described.


Recommended Dilution

ELISA: 1:10000

WB: 1:500 - 1:2000

Not yet tested in other applications.

Images



Western Blot analysis using GTBP Monoclonal antibody against truncated MSH6 recombinant protein.

Storage

-20°C for one year