

KLHL25 Monoclonal Antibody

Description

Product type	Primary Antibody
Code	BT-MCA0818
Host	Mouse
Isotype	IgG
Size	50ul, 100ul
Immunogen	Purified recombinant fragment of human KLHL25 expressed in E. Coli.
Mol wt	N/A
Species reactivity	Human
Clonality	Monoclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	Ectoderm-neural cortex protein 2
Synonyms	KLHL25; ENC2; Ectoderm-neural cortex protein 2; ENC-2; Kelch-like protein 25

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

N/A


Recommended Dilution

ELISA: 1:10000

WB: 1:500 - 1:2000

Not yet tested in other applications.

Images



Western Blot analysis using KLHL25 Monoclonal antibody against KLHL25-hIgGFc transfected HEK293 cell.

Storage

-20°C for one year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com