

MUM1 Monoclonal Antibody

Description

Product type	Primary Antibody
Code	BT-MCA0905
Host	Mouse
Isotype	IgG
Size	50ul, 100ul
Immunogen	Purified recombinant fragment of human MUM1 expressed in E. Coli.
Mol wt	N/A
Species reactivity	Human
Clonality	Monoclonal
Recommended application	WB, IHC-p, IF, ELISA
Concentration	1 mg/ml
Full name	PWWP domain-containing protein MUM1
Synonyms	MUM1; EXPAND1; PWWP domain-containing protein MUM1; Mutated melanoma-associated antigen 1; MUM-1; Protein expandere

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,similarity:Belongs to the MUM1 family.,

Recommended Dilution


ELISA: 1:10000

IHC: 1:200 - 1:1000


WB: 1:500 - 1:2000

Not yet tested in other applications.

Images



Western Blot analysis using MUM1 Monoclonal antibody against recombinant fragment of human MUM1.



Immunohistochemistry analysis of paraffin-embedded human melanoma (A), non-Hodgkin's lymphoma (B), with DAB staining using MUM1 Monoclonal antibody.

Storage

-20°C for one year