

PEG10 Monoclonal Antibody

Description

Product type	Primary Antibody
Code	BT-MCA1017
Host	Mouse
Isotype	IgG
Size	50ul, 100ul
Immunogen	Purified recombinant human PEG10 protein fragments expressed in E.coli.
Mol wt	N/A
Species reactivity	Human,Mouse,Pig,Rabbit
Clonality	Monoclonal
Recommended application	WB
Concentration	1 mg/ml
Full name	Retrotransposon-derived protein PEG10
Synonyms	PEG10; EDR; KIAA1051; MAR2; MART2; MEF3L1; RGAG3; Retrotransposon-derived protein PEG10; Embryonal carcinoma differentiation-regulated protein; Mammalian retrotransposon-derived protein 2; Myelin expression factor 3-like protein 1; MEF3-lik

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background


This is a paternally expressed imprinted gene that is thought to have been derived from the Ty3/Gypsy family of retrotransposons. It contains two overlapping open reading frames, RF1 and RF2, and expresses two proteins: a shorter, gag-like protein (with a CCHC-type zinc finger domain) from RF1 and a longer, gag/pol-like fusion protein (with an additional aspartic protease motif) from RF1/RF2 by -1 translational frameshifting (-1 FS). While -1 FS has been observed in RNA viruses and transposons in both prokaryotes and eukaryotes, this gene represents the first example of -1 FS in a eukaryotic cellular gene. This gene is highly conserved across mammalian species and retains the heptanucleotide (GGGAAAC) and pseudoknot elements required for -1 FS. It is expressed in adult and embryonic tissues (most notably in placenta) and reported to have a role in cell proliferation, differentiation and apoptosis.

Recommended Dilution

WB: 1:1000 - 1:2000

Not yet tested in other applications.

Images



Western Blot analysis using PEG10 Monoclonal antibody against HeLa cell lysate.

Storage

-20°C for one year