

Tyro3 Monoclonal Antibody

Description

Product type	Primary Antibody
Code	BT-MCA1260
Host	Mouse
Isotype	IgG
Size	50ul, 100ul
Immunogen	Purified recombinant fragment of Tyro3 expressed in E. Coli.
Mol wt	N/A
Species reactivity	Human
Clonality	Monoclonal
Recommended application	WB, IF, ICC, ELISA
Concentration	1 mg/ml
Full name	Tyrosine-protein kinase receptor TYRO3
Synonyms	TYRO3; BYK; DTK; RSE; SKY; Tyrosine-protein kinase receptor TYRO3; Tyrosine-protein kinase DTK; Tyrosine-protein kinase RSE; Tyrosine-protein kinase SKY; Tyrosine-protein kinase byk

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

The gene is part of a 3-member transmembrane receptor kinase receptor family with a processed pseudogene distal on chromosome 15. The encoded protein is activated by the products of the growth arrest-specific gene 6 and protein S genes and is involved in controlling cell survival and proliferation, spermatogenesis, immunoregulation and phagocytosis. The encoded protein has also been identified as a cell entry factor for Ebola and Marburg viruses.

Recommended Dilution


ELISA: 1:10000

IF: 1:200 - 1:1000

WB: 1:500 - 1:2000

Not yet tested in other applications.

Images



Western Blot analysis using Tyro3 Monoclonal antibody against truncated Tyro3 recombinant protein.

Storage

-20°C for one year