

## MAP2K4 Monoclonal Antibody

### Description

Product type	Antibody
Code	BT-MCA2366
Host	Mouse
Isotype	Mouse IgM
Size	100µL, 50µL
Immunogen	Purified recombinant fragment of MAP2K4 expressed in E. Coli.
Mol wt	N/A
Species reactivity	Human
Clonality	Monoclonal
Recommended application	IHC
Concentration	N/A
Full name	N/A
Synonyms	MAP2K4

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

### Background

N/A


### Recommended Dilution

IHC-p: 1:200 - 1:1000

ELISA: 1:10000

Not yet tested in other applications.

### Images



Immunohistochemical analysis of paraffin-embedded human thymoma tissue (left) and spleen tissue (right), showing cytoplasmic localization using MAP2K4 mouse mAb with DAB staining.

### Storage

Store at 4°C short term. Aliquot and store at -20°C long term.

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: [save@bt-laboratory.com](mailto:save@bt-laboratory.com) | [www.bt-laboratory.com](http://www.bt-laboratory.com)