

GSTM1 Monoclonal Antibody

Description

Product type	Antibody
Code	BT-MCA3648
Host	Mouse
Isotype	Mouse IgG1
Size	100µL, 50µL
Immunogen	Purified recombinant fragment of human GSTM1 expressed in E. Coli.
Mol wt	26kDa
Species reactivity	Human,Monkey
Clonality	Monoclonal
Recommended application	WB,IHC,FCM
Concentration	N/A
Full name	N/A
Synonyms	MU;H-B;GST1;GTH4;GTM1;MU-1;GSTM1-1;MGC26563;GSTM1a-1a;GSTM1b-1b

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

Cytosolic and membrane-bound forms of glutathione S-transferase are encoded by two distinct supergene families. At present, eight distinct classes of the soluble cytoplasmic mammalian glutathione S-transferases have been identified: alpha, kappa, mu, omega, pi, sigma, theta and zeta. This gene encodes a glutathione S-transferase that belongs to the mu class. The mu class of enzymes functions in the detoxification of electrophilic compounds, including carcinogens, therapeutic drugs, environmental toxins and products of oxidative stress, by conjugation with glutathione. The genes encoding the mu class of enzymes are organized in a gene cluster on chromosome 1p13.3 and are known to be highly polymorphic. These genetic variations can change an individual's susceptibility to carcinogens and toxins as well as affect the toxicity and efficacy of certain drugs. Null mutations of this class mu gene have been linked with an increase in a number of cancers, likely due to an increased susceptibility to environmental toxins and carcinogens. Multiple protein isoforms are encoded by transcript variants of this gene.

Recommended Dilution

WB: 1:500 - 1:2000


IHC-p: 1:200 - 1:1000

FCM: 1:200 - 1:400


ELISA: 1:10000

Not yet tested in other applications.


Images




Western blot analysis using GSTM1 mAb against human GSTM1 (AA: 23-181) recombinant protein.
(Expected MW is 25.7 kDa)




Western blot analysis using GSTM1 mouse mAb against Cos7 (1), MCF-7 (2), Jurkat (3), HeLa (4), HL7702 (5) and HepG2 (6) cell lysate.




Immunohistochemical analysis of paraffin-embedded rectum cancer tissues using GSTM1 mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded liver cancer tissues using GSTM1 mouse mAb with DAB staining.



Flow cytometric analysis of Jurkat cells using GSTM1 mouse mAb (green) and negative control (red).



Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);

Storage

Store at 4°C short term. Aliquot and store at -20°C long term.

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com