

CD30 Monoclonal Antibody

Description

| | |
|--------------------------------|---|
| Product type | Antibody |
| Code | BT-MCA3772 |
| Host | Mouse |
| Isotype | Mouse IgG1 |
| Size | 100µL, 50µL |
| Immunogen | Purified recombinant fragment of human CD30 (AA: 416-595) expressed in E. Coli. |
| Mol wt | 63.7kDa |
| Species reactivity | Human,Mouse |
| Clonality | Monoclonal |
| Recommended application | WB,ICC,FCM |
| Concentration | N/A |
| Full name | N/A |
| Synonyms | Ki-1;D1S166E,TNFRSF8 |

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

The protein encoded by this gene is a member of the TNF-receptor superfamily. This receptor is expressed by activated, but not by resting, T and B cells. TRAF2 and TRAF5 can interact with this receptor, and mediate the signal transduction that leads to the activation of NF-kappaB. This receptor is a positive regulator of apoptosis, and also has been shown to limit the proliferative potential of autoreactive CD8 effector T cells and protect the body against autoimmunity. Two alternatively spliced transcript variants of this gene encoding distinct isoforms have been reported.

Recommended Dilution

WB: 1:500 - 1:2000


ICC: 1:50 - 1:200

FCM: 1:200 - 1:400


ELISA: 1:10000

Not yet tested in other applications.


Images




Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)




Western blot analysis using CD30 mAb against human CD30 (AA: 416-595) recombinant protein. (Expected MW is 22.4 kDa)




Western blot analysis using CD30 mAb against HEK293-6e (1) and CD30 (AA: 416-595)-hIgGFc transfected HEK293-6e (2) cell lysate.



Western blot analysis using CD30 mouse mAb against HeLa (1), MCF-7 (2), HCT116 (3), and PC-3 (4) cell lysate.



Immunofluorescence analysis of HeLa cells using CD30 mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye. Red: Actin filaments have been labeled with Alexa Fluor- 555 phalloidin. Secondary antibody from Fisher (Cat#: 35503)



Flow cytometric analysis of Jurkat cells using CD30 mouse mAb (green) and negative control (red).

Storage

Store at 4°C short term. Aliquot and store at -20°C long term.

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com