

VISTA Monoclonal Antibody

Description

Product type	Antibody
Code	BT-MCA4087
Host	Mouse
Isotype	Mouse IgG1
Size	100µL, 50µL
Immunogen	Purified recombinant fragment of human VISTA (AA: extra 33-194) expressed in HEK293 cells.
Mol wt	34kDa
Species reactivity	Human
Clonality	Monoclonal
Recommended application	FCM
Concentration	N/A
Full name	N/A
Synonyms	VSIR;B7H5;GI24;B7-H5;PD-1H;SISP1;PP2135;C10orf54;DD1alpha

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

N/A

Recommended Dilution


WB: 1:500 - 1:2000

FCM: 1:200 - 1:400


ELISA: 1:10000

Not yet tested in other applications.


Images



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Western blot analysis using VISTA mAb against human VISTA (AA: extra 33-194) recombinant protein. (Expected MW is 48 kDa)



Flow cytometric analysis of Jurkat cells using VISTA mouse mAb (green) and negative control (red).

Storage

Store at 4°C short term. Aliquot and store at -20°C long term.

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com