

Phospho-APLF (S116) Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-PHS01078
Host	Rabbit
Isotype	IgG
Size	20ul, 50ul, 100ul
Immunogen	The antiserum was produced against synthesized peptide derived from human APLF around the phosphorylation site of Ser116. AA range:82-131
Mol wt	56956
Species reactivity	Human, mouse
Clonality	Polyclonal
Recommended application	IHC-p, ELISA
Concentration	1 mg/ml
Full name	Phospho-APLF (S116) Antibody
Synonyms	APLF; C2orf13; PALF; XIP1; Aprataxin and PNK-like factor; Apurinic-apyrimidinic endonuclease APLF; PNK and APTX-like FHA domain-containing protein; XRCC1-interacting protein 1

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

C2ORF13 is a component of the cellular response to chromosomal DNA single- and double-strand breaks (Iles et al. 2007 [PubMed 17353262]).

Recommended Dilution

IHC: 1: 100 - 1: 300

ELISA: 1: 5000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for one year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com