

Phospho-53BP1 (S25) Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-PHS01140
Host	Rabbit
Isotype	IgG
Size	20ul, 50ul, 100ul
Immunogen	The antiserum was produced against synthesized peptide derived from human 53BP1 around the phosphorylation site of Ser25. AA range:10-59
Mol wt	213574
Species reactivity	Human, mouse, rat
Clonality	Polyclonal
Recommended application	IHC-p, IF, ELISA
Concentration	1 mg/ml
Full name	Phospho-53BP1 (S25) Antibody
Synonyms	TP53BP1; Tumor suppressor p53-binding protein 1; 53BP1; p53-binding protein 1; p53BP1

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

TP53BP1 (Tumor Protein P53 Binding Protein 1) is a Protein Coding gene. Diseases associated with TP53BP1 include riddle syndrome and ataxia-telangiectasia. Among its related pathways are DNA Damage and ATM Pathway. GO annotations related to this gene include sequence-specific DNA binding and p53 binding.

Recommended Dilution

IHC: 1: 100 - 1: 300

IF: 1: 200 - 1: 1000

ELISA: 1: 10000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for one year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com