



Optimize Your Research

# Recombinant Ovine Interferon-tau protein

## Code CD00156

**Storage:** This lyophilized preparation is stable at 2-8°C, but should be kept at -20°C for long term storage, preferably desiccated. Upon reconstitution, the preparation is stable for up to one week at 2-8°C. For maximal stability, apportion the reconstituted preparation into working aliquots and store at -20°C to -70°C. Avoid repeated freeze/thaw cycles.

### Intended Use

This product is for research use only, not for use in diagnosis procedures. It is highly recommended to read this instruction entirely before the use.

### Source

Pichia Pastoris

### Molecular Weight

Approximately 19.9 kDa, a single glycosylated polypeptide chain containing 172 amino acids.

### Purity

> 97 % by SDS-PAGE and HPLC analyses.

### Biological Activity

Fully biologically active when compared to IFN-alpha. The specific activity determined by a viral resistance assay is no less than  $1.0 \times 10^7$  IU/mg.

### Physical Appearance

Sterile Filtered White lyophilized (freeze-dried) powder.

### Formulation

Lyophilized from a 0.2 µm filtered concentrated solution in PBS, pH 7.4.

### Sequence

CYLSRKMLLD ARENLKLLDR MNRLSPHSCL QDRKDFGLPQ EMVEGDQLQK DQAFPVLYEM LQQSFNLFYT EHSSAAWDTT  
LLEQLCTGLQ QQLDHLDTCR GQVMGEEDSE LGNMDPIVTV KKYFQGIYDY LQEKGYSDCA WEIVRVEMMR ALTVSTTLQK  
RLTKMGGDLN SP

### Endotoxin

Less than 0.1 EU/µg of rOvIFN-τ as determined by LAL method.

### Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Reconstitute in sterile distilled water or aqueous buffer containing 0.1 % BSA to a concentration of 0.1-1.0 mg/mL. Stock solutions should be apportioned into working aliquots and stored at ≤ -20°C. Further dilutions should be made in appropriate buffered solutions.

If you have any question on order please contact us via: [order@bt-laboratory.com](mailto:order@bt-laboratory.com); technical assistance please contact us via: [support@bt-laboratory.com](mailto:support@bt-laboratory.com) More product visit [www.bt-laboratory.com](http://www.bt-laboratory.com)

5/F 2 Bldg., 501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China  
Tel: 86 21 31007137 | E-mail: [save@bt-laboratory.com](mailto:save@bt-laboratory.com) | [www.bt-laboratory.com](http://www.bt-laboratory.com)