



Optimize Your Research

# Recombinant Human Interleukin-36 alpha, 158a.a. protein

## Code CD01505

**Storage:** This lyophilized preparation is stable at 2-8°C, but should be kept at -20°C for long term storage, preferably desiccated. Upon reconstitution, the preparation is stable for up to one week at 2-8°C. For maximal stability, apportion the reconstituted preparation into working aliquots and store at -20°C to -70°C. Avoid repeated freeze/thaw cycles.

### Intended Use

This product is for research use only, not for use in diagnosis procedures. It is highly recommended to read this instruction entirely before the use.

### Source

Escherichia coli.

### Molecular Weight

Approximately 17.7 kDa, a single non-glycosylated polypeptide chain containing 158 amino acids.

### Purity

> 95 % by SDS-PAGE and HPLC analyses.

### Biological Activity

Fully biologically active when compared to standard. The specific activity determined by its ability in a functional ELISA. Immobilized rHuIL-36 $\alpha$  at 1  $\mu$ g/mL can bind recombinant human IL-1 Rrp2 Fc Chimera with a range of 0.15-5  $\mu$ g/mL.

### Physical Appearance

Sterile Filtered White lyophilized (freeze-dried) powder.

### Formulation

Lyophilized from a 0.2  $\mu$ m filtered concentrated solution in 2  $\times$  PBS, pH 7.4.

### Sequence

MEKALKIDTP QQGSIQDINH RVWVLQDQTL IAVPRKDRMS PVTIALISCR HVETLEKDRG NPIYLGLNGL NLCLMCAKVG DQPTLQLKEK DIMDLYNQPE PVKSFLFYHS QSGRNSTFES VAFPGWFIIV SSEGGCPLIL TQELGKANTT DFGLTMLF

### Endotoxin

Less than 1 EU/ $\mu$ g of rHuIL-36 $\alpha$ , 158a.a. as determined by LAL method.

### Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Reconstitute in sterile distilled water or aqueous buffer containing 0.1 % BSA to a concentration of 0.1-1.0 mg/ml. Stock solutions should be apportioned into working aliquots and stored at  $\leq$  -20 C. Further dilutions should be made in appropriate buffered solutions.

If you have any question on order please contact us via: [order@bt-laboratory.com](mailto:order@bt-laboratory.com); technical assistance please contact us via: [support@bt-laboratory.com](mailto:support@bt-laboratory.com) More product visit [www.bt-laboratory.com](http://www.bt-laboratory.com)