



Optimize Your Research

Recombinant Human Insulin-like Growth Factor-Binding Protein 7 protein

Code CD01586

Storage: This lyophilized preparation is stable at 2-8°C, but should be kept at -20°C for long term storage, preferably desiccated. Upon reconstitution, the preparation is stable for up to one week at 2-8°C. For maximal stability, apportion the reconstituted preparation into working aliquots and store at -20°C to -70°C. Avoid repeated freeze/thaw cycles.

Intended Use

This product is for research use only, not for use in diagnosis procedures. It is highly recommended to read this instruction entirely before the use.

Source

Escherichia coli.

Molecular Weight

Approximately 26.4 kDa, a single non-glycosylated polypeptide chain containing 256 amino acids.

Purity

> 95 % by SDS-PAGE and HPLC analyses.

Biological Activity

Testing in Progress.

Physical Appearance

Sterile Filtered White lyophilized (freeze-dried) powder.

Formulation

Lyophilized from a 0.2 µm filtered concentrated solution in 20 mM Tris-HCl, pH 8.6, 150 mM NaCl.

Sequence

SSSDTCGPCE PASCPLPPL GLLGETRDA CGCCPMCARG EGEPGGGGA GRGYCAPGME CVKSRKRRKG KAGAAAGGPG
VSGVCVCKSR YPVCSDGTT YPSGCQLRAA SQRAESRGEK AITQVSKGTC EQGPSIVTTP KDIWNVTGAQ VYLSCEVIGI PTPVLIWNKV
KRGHYGVQRT ELLPGDRDNL AIQTRGGPEK HEVTGWVLVS PLSKEDAGEY ECHASNSQGQ AS

Endotoxin

Less than 0.1 EU/µg of rHuIGF-BP7 as determined by LAL method.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Reconstitute in sterile distilled water or aqueous buffer containing 0.1 % BSA to a concentration of 0.1-1.0 mg/mL. Stock solutions should be apportioned into working aliquots and stored at ≤ -20 °C. Further dilutions should be made in appropriate buffered solutions.

If you have any question on order please contact us via: order@bt-laboratory.com; technical assistance please contact us via: support@bt-laboratory.com More product visit www.bt-laboratory.com