



Optimize Your Research

Recombinant Human Fibroblast Growth Factor-19 protein

Code CD01806

Storage: This lyophilized preparation is stable at 2-8°C, but should be kept at -20°C for long term storage, preferably desiccated. Upon reconstitution, the preparation is stable for up to one week at 2-8°C. For maximal stability, apportion the reconstituted preparation into working aliquots and store at -20°C to -70°C. Avoid repeated freeze/thaw cycles.

Intended Use

This product is for research use only, not for use in diagnosis procedures. It is highly recommended to read this instruction entirely before the use.

Source

Escherichia coli.

Molecular Weight

Approximately 21.8 kDa, a single non-glycosylated polypeptide chain containing 195 amino acids.

Purity

> 95 % by SDS-PAGE and HPLC analyses.

Biological Activity

Assay #1: Fully biologically active when compared to standard. The biological activity is measured by its binding ability in a functional ELISA. Immobilized rHuFGF R4 at 5 µg/ml can bind rHuFGF-19 with a linear range of 3-200 ng/ml. Assay #2: Fully biologically active when compared to standard. The ED50 as determined by a cell proliferation assay using Mouse Balb/c 3T3 cells is less than 150 ng/ml, corresponding to a specific activity of > 6.7 × 10³ IU/mg.

Physical Appearance

Sterile Filtered White lyophilized (freeze-dried) powder.

Formulation

Lyophilized from a 0.2 µm filtered concentrated solution in PBS, pH 7.4.

Sequence

MRPLAFSDAG PHVHYGWGDP IRLRHLYTSG PHGLSSCFLLR IRADGVVDCA RGQSAHSLE IKAVALRTVA IKGVHSVRYL
CMGADGKMQG LLQYSEEDCA FEEIIRPDGY NVYRSEKHRL PVSLSSAKQR QLYKNRGFLP LSHFLPMLPM VPEEPEDLRG
HLESDMFSSP LETDSMDPFG LVTGLEAVRS PSFEK

Endotoxin

Less than 1 EU/µg of rHuFGF-19 as determined by LAL method.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Reconstitute in sterile distilled water or aqueous buffer containing 0.1 % BSA to a concentration of 0.1-1.0 mg/mL. Stock solutions should be apportioned into working aliquots and stored at ≤ -20 °C. Further dilutions should be made in appropriate buffered solutions.

If you have any question on order please contact us via: order@bt-laboratory.com; technical assistance please contact us via: support@bt-laboratory.com More product visit www.bt-laboratory.com