

Optimize Your Research

Recombinant Human Cu/Zn Superoxide Dismutase, His protein

Code CD01967

Storage: This lyophilized preparation is stable at 2-8°C, but should be kept at -20°C for long term storage, preferably desiccated. Upon reconstitution, the preparation is stable for up to one week at 2-8°C. For maximal stability, apportion the reconstituted preparation into working aliquots and store at -20°C to -70°C. Avoid repeated freeze/thaw cycles.

Intended Use

This product is for research use only, not for use in diagnosis procedures. It is highly recommended to read this instruction entirely before the use.

Source

Escherichia coli.

Molecular Weight

Approximately 39.9 kDa, a homodimer, non-glycosylated polypeptide chain containing 2 × 189 amino acids with Met, Gly and 10 × His at N-terminus.

Purity

> 95 % by SDS-PAGE and HPLC analyses.

Biological Activity

Fully biologically active when compared to standard. The potency per mg was determined by Pyrogalllic Acid method and was found to be more than 1.0 × 10⁴ IU/mg.

Physical Appearance

Sterile Filtered White lyophilized (freeze-dried) powder.

Formulation

Lyophilized from a 0.2 µm filtered concentrated solution in PBS, pH 7.4.

Sequence

MGHHHHHHHHH HHSSGHIEGR HMTYARAAAR QARALEATKA VCVLKGDGPV QGIINFEQKE SNGPVKWWGS IKGLTEGLHG
FHVHEFGDNT AGCTSAGPHF NPLSRKHGGP KDEERHVGDL GNVTADKDG V ADVSIEDSVI SLSGDHCIIG RTLVVHEKAD
DLGKGGNEES TKTGNAGSRL ACGVIGIAQ

Endotoxin

Less than 1 EU/µg of rHuCu/Zn SOD, His as determined by LAL method.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Reconstitute in sterile distilled water or aqueous buffer containing 0.1 % BSA to a concentration of 0.1-1.0 mg/mL. Stock solutions should be apportioned into working aliquots and stored at ≤ -20 °C. Further dilutions should be made in appropriate buffered solutions.

If you have any question on order please contact us via: order@bt-laboratory.com; technical assistance please contact us via: support@bt-laboratory.com More product visit www.bt-laboratory.com